EN 7 I.U.M.C. Endometrial Sampler

Instructions for Use

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Komora 8 ml
Câmara de 8 ml
8 ml kammare

15. For additional testing or teaching slides

Pro další testovací nebo výuková sklíčka
Til ekstra test- eller undervisningsobjektglas
Für weitere Tests oder Lehr-Objektgläser
Για πρόσθετες δοκιμαστικές ή εκπαιδευτικές αντικειμενοφόρες πλάκες
Para pruebas adicionales o para preparar portaobjetos para enseñanza
Pour analyses supplémentaires ou lames pédagogiques
További tesztelési vagy oktatási célú üveglapokhoz
Vetrini per ulteriori analisi o a scopo didattico
Voor extra test- of onderwijsobjectglazen
For objektglass til ytterligere testing eller undervisning
Na dodatkowe badania lub preparaty szkoleniowe
Para testes adicionais ou lâminas para ensino
För ytterligare test eller objektglas för undervisning

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I.U.M.C. ENDOMETRIAL SAMPLER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Set Contains (*set components may vary*)

- Stainless steel shaft, 20 cm long with nylon brush and plastic handle with reference mark
- 9.0 French vinyl sheath, 16 cm long

INTENDED USE

Used to obtain endometrial cell samples.

NOTE: The most reliable method for the detection of endometrial carcinoma is a formal dilation of the cervix, endometrial curettage and exploration of the uterine cavity under general anesthesia. However, there are situations when a simpler sampling technique is adequate, and the healthcare provider is willing to accept the limitations of the sampling technique. This endometrial sampler is intended for use in various endometrial cytological cell sampling procedures including early detection of occult endometrial carcinoma among asymptomatic patients for whom curettage is not indicated, e.g., screening of patients at risk of endometrial carcinoma, follow-up of patients receiving estrogen replacement therapy or Tamoxifen, follow-up of treated and untreated patients with endometrial hyperplasia, as well as evaluation of postmenopausal women whose cervico-vaginal smears show a high estrogen effect and endometrial cells.

CONTRAINDICATIONS

This device should not be used during pregnancy or following a recent uterine perforation. In patients with a history of Cesarean section, use this device with caution. The Tao Brush™ is not intended to be used if the cervix is infected or inflamed.

WARNINGS

None known

PRECAUTIONS

- In patients with postmenopausal bleeding who are at known increased risk of carcinoma, and when an inadequate tissue sample has been removed by a biopsy device, consider 1) obtaining additional transcervical cytologic samples by an endometrial sampler in a physician's office, or, 2) dilation and curettage under general anesthesia in a hospital.
- The patient should experience only mild discomfort and little or no pain during the procedure. If the patient experiences significant pain, the endometrial sampler should be withdrawn immediately.
- Insertion of a vaginal speculum and dilation of the cervix is unnecessary; however, if the brush and sheath do not pass into the endometrial cavity with gentle pressure, it is possible that there is cervical stenosis, and dilation and curettage under general anesthesia should be considered.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

Using The Tao Brush™

1. Position screw-cap test tube containing 8 ml of CytoRich® Brush Cytology Preservative (AutoCyte, Inc., Elon College, NC) in a test-tube rack at the site of the procedure.
2. Place patient in lithotomy position.
3. Retract the brush sampler completely into the outer sheath. (**Figure A**)

4. Gently insert the device to the level of the fundus. **(Figure B)**
 5. Pull back the outer sheath all the way to the handle. Amply rotate the brush sampler. **(Figure C)** Two methods are suggested: 1) Rotate brush sampler in a clockwise manner until reference mark on the handle indicates completion of a 360° turn, then rotate counterclockwise the (opposite direction) until the reference mark on the handle indicates completion of a 360° turn; 2) Rotate the brush sampler in only one direction by completing 4 or 5 360° rotations. **NOTE:** Reference mark on handle indicates completion of a 360° rotation.
 6. In order to trap endometrial material *in situ*, push the outer sheath over the brush to the tip and remove the device. **(Figure D)** The normal endometrial cavity is in a collapsed state, so the brush will have direct contact with the entire endometrial surface.
 7. Immediately immerse the device into 8 ml of CytoRich® Brush Cytology Preservative.
 8. Retract the sheath to expose the brush to preservative solution.
 9. Hold the sheath firmly and move the brush in and out of the sheath to clean it of adherent cells and tissue. **(Figure E)**
- NOTE:** Collections are stable in preservative for periods of up to several weeks.
10. Remove the brush assembly from the test-tube, replace the screw cap, and submit the tube to the laboratory for processing.

To Obtain Uncontaminated Endometrial Cultures

1. After insertion of a sterile, nonlubricated vaginal speculum, swab the ectocervix and the endocervical canal with povidone iodine solution. **NOTE:** Insert the swab about 1.5 cm into the endocervical canal to ensure adequate swabbing of the endocervix with the povidone.
2. Insert the brush into the endometrial cavity following steps 3-6 from the INSTRUCTIONS FOR USE section preceding these instructions. The reference mark on the handle indicates completion of a 360° turn.
3. Remove sampler.
4. Wipe the rounded tip of the brush with 95% alcohol gauze.
5. Pull back the sheath. Prepare morphologic evaluation (if required) by preparing a direct smear on a sterile glass slide and spray-fix immediately.
6. For culture studies, place the brush into sterile Stuarts Transportation Medium and agitate for 5 seconds.

Lab Instructions For Sample Processing

NOTE: This process will yield several sample slides with an 18.5 mm diameter circle of cell material. (Generally, only two slides are prepared for diagnostic purposes.) Experience shows that because cells and their aggregates are homogeneously suspended in the dispersal solution, only representative slides are needed for diagnosis.

NOTE: Pre-coat glass slides used in processing with poly-L-lysine (CytoRich® Slide Coat Reagent; AutoCyte Inc., Elon College, NC). Pre-soaking the glass slides in a weak (e.g., 0.01 Normal) ammonium hydroxide solution for about 10 minutes (rinse in deionized water) before coating them with poly-L-lysine enhances cell material adherence to the slide and facilitates slide storage.

1. Vortex the CytoRich® brush cytology tube to mix well.
2. Place in a test rack for 1 minute to allow tissue fragments to settle to the bottom of the tube.
3. Using a polyethylene transfer pipette, transfer the upper 6 ml of sample to a clean 50 ml centrifuge tube with cap.
4. Transfer the microbiopsies remaining in the residual suspension bottom 2 ml to formalin fixative for cell block processing, using Histoscreen® (Perk® Scientific Inc., Yeadon, PA) and submit to histology lab.
5. Add 30 ml of CytoRich® Yellow to the 6 ml cytology suspension in the 50 ml centrifuge tube.

6. Allow the suspension to equilibrate for 10 minutes.
7. Use the Hettich Universal Cytocentrifuge (AutoCyte Inc., Elon College, NC) fitted with 8 ml sample chambers to prepare slides.
8. Depending on sample cellularity, add 3 to 5 ml of the suspension from step 5 to each of the two chambers.
9. Set the time switch to 15 minutes.
10. Set the speed regulator to 30%.
11. After cytocentrifugation, aspirate supernate off the monolayer preparation and disassemble sample chambers.
12. Air dry the slides for 15-20 minutes.
13. Post fix the slides in 95% ethanol for 10 minutes.
14. Wash in warm running water for 1-2 minutes.
15. Stain using the Papanicolaou stain or its variants.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCE

J. A. Maksem: "Ciliated Cell Adenocarcinoma of the Endometrium Diagnosed by Endometrial Brush Cytology and Confirmed by Hysterectomy: A Case Report Detailing a Highly Efficient Cytology Collection and Processing Technique," Diagnostic Cytopathology, 16 (1997), 78-82.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.