EN 3	FluoroSet™ Radiographic Tubal Assessment Set Instruction for Use
CS 4	FluoroSet™ Souprava pro rentgenové vyšetření vejcovodů Návod k použití
DA 6	FluoroSet™ Radiografisk tubavurderingssæt Brugsanvisning
DE 7	FluoroSet™ Set zur radiographischen Beurteilung der Eileiter Gebrauchsanweisung
EL 9	FluoroSet™ Σετ ακτινογραφικής αξιολόγησης σαλπίνγων Οδηγίες χρήσης
ES 11	FluoroSet™ Equipo para la valoración radiográfica de la función tubárica Instrucciones de uso
FR 12	FluoroSet™ Set pour évaluation radiographique tubaire Mode d'emploi
HU 14	FluoroSet™ Radiográfiás petevezeték-vizsgáló készlet Használati utasítás
IT 16	FluoroSet™ Set di valutazione radiografica delle tube Istruzioni per l'uso
NL 17	FluoroSet™ Set voor radiografisch tubair onderzoek Gebruiksaanwijzing
NO 19	FluoroSet™ Radiografisk tubarvurderingssett Bruksanvisning
PL 21	FluoroSet™ Zestaw do radiograficznej oceny jajowodów Instrukcja użycia
PT 22	FluoroSet™ Conjunto para avaliação radiográfica das trompas de Falópio Instruções de utilização
SV 24	FluoroSet™ Set för äggledarbedömning med röntgen Bruksanvisning

FLUOROSSET™

RADIOGRAPHIC TUBAL ASSESSMENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Set Contains (*set components may vary*)
- Intrauterine Access Balloon Catheter
- Selective Salpingography Catheter with Beacon® Tip
- Roadrunner PC Wireguide with SlipCoat Hydrophilic Coating

INTENDED USE

Used for instillation of contrast media into the uterine cavity for radiographic evaluation of the uterine cavity and for injection of appropriate contrast media into the fallopian tubes for evaluation of tubal patency

CONTRAINDICATIONS

Do not use in a patient with an active vaginal or intrauterine infection, a sexually transmitted disease, a recent uterine perforation, a recent pregnancy (or is currently pregnant), or if the patient currently has an intrauterine device.

WARNINGS

- Due to thin wall construction of the Selective Salpingography Catheter, extreme care must be exercised during manipulation and withdrawal to prevent pulling the catheter apart. Reshaping of the tip is not recommended. Damage can result when exposed to heat.
- The FluoroSet™ Radiographic Tubal Assessment Set is not intended for complete catheterization of the fallopian tube or tubal recanalization.
- This catheter system is only designed to access the fallopian tube from the uterine cavity. It should only be advanced as far as necessary to achieve a meaningful diagnostic study. Possible adverse effects of more distal advancement have not yet been determined. The catheter system is not intended for tubal recanalization, and there is no data available, to date, showing any clinical benefit when used this way.

PRECAUTIONS

- Do not exceed the recommended balloon volume. Over-inflation of this device can cause the balloon to rupture. Inflate the balloon with only sterile saline. Do not inflate the balloon with air.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Advance Intrauterine Access Balloon Catheter transcervically until the balloon is positioned in the lower portion of the cervix or uterus. Inflate the balloon to maximum inflation volume as stated on the product label or on the inflation check valve.
2. Attach a Luer Lock syringe to the sidearm assembly and stopcock to facilitate instillation of 5-10 mL of an appropriate contrast media into uterus for basic hystero-graphy. This volume of fluid is to outline the basic anatomy and provide guidance in preparation for Selective Salpingography (SSG) Catheter placement.
3. Insert Selective Salpingography (SSG) Catheter through the working channel of the Intrauterine Access Balloon catheter until the ink mark on the SSG Catheter enters the Check-Flo Adapter located on the

proximal end of the Intrauterine Access Balloon. While scanning with fluoroscopy, rotate catheter to address the appropriate fallopian tube. The distal tip of the SSG Catheter is radiopaque and will be visible under fluoroscopy. An appropriate contrast media may now be injected through the SSG Catheter to assess tubal patency.

4. If access to the tubes is hindered or difficult, the Roadrunner PC Wire Guide may be used to facilitate tubal access.

INSTRUCTIONS FOR USING THE ROADRUNNER PC WIRE GUIDE

ACTIVATING HYDROPHILIC COATING

Roadrunner PC Wire Guides are very slippery when wet. Always maintain control of the wire guide with the enclosed Olcott manipulation device when maneuvering it through any access device.

1. Prior to using the wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder.
2. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating.
 - For optimal performance, after exposure to ambient environment, it is necessary to rehydrate the wire guide prior to instrument exchange. Hydrophilic coating is not permanent. If after extended use the coated wire guide does not perform smoothly, replace the wire guide with a new hydrophilic wire guide.
3. Introduce the Roadrunner PC Wire Guide through the Selective Salpingography (SSG) Catheter and into the uterine cavity. The wire guide should be visible under radiographic guidance. Use the wire guide to facilitate fallopian tube access at the level of the fundal ostia.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.