- EN**  
3  
Russell Gastrostomy Tray with Cope  
Gastrointestinal Suture Anchor Set  
Instructions for Use
- DA**  
4  
Russell gastrostomibakke med Cope  
gastrointestinale suturankersæt  
Brugsanvisning
- DE**  
6  
Russell Gastrostomie-Komplettset mit Cope  
gastrointestinalem Nahtanker-Set  
Gebrauchsanweisung
- EL**  
8  
Δίσκος γαστροστομίας Russell με σετ  
γαστρεντερικού αγκίστρου ράμματος Cope  
Οδηγίες χρήσης
- ES**  
10  
Bandeja de gastrostomía Russell con equipo  
de anclajes de sutura gastrointestinal Cope  
Instrucciones de uso
- FR**  
12  
Plateau de gastrostomie Russell avec set  
d'ancre de suture gastro-intestinale Cope  
Mode d'emploi
- IT**  
14  
Set completo per gastrostomia Russell con set  
di suture con ancora gastrointestinale Cope  
Istruzioni per l'uso
- NL**  
16  
Russell-gastrostomiepakket met Cope  
gastro-intestinale hechtankerset  
Gebruiksaanwijzing
- PT**  
18  
Tabuleiro de gastrostomia Russell com  
conjunto de ancoragem de sutura  
gastrointestinal Cope  
Instruções de utilização
- SV**  
20  
Russell gastrostomibricka med Cope  
gastrointestinale suturankersæt  
Bruksanvisning

## RUSSELL GASTROSTOMY TRAY WITH COPE GASTROINTESTINAL SUTURE ANCHOR SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Russell Gastrostomy Tray consists of a silicone gastrostomy catheter, Peel-Away® Introducer, and appropriately sized introducer needle, wire guide, and dilators. The tray also contains a Cope Gastrointestinal Suture Anchor Set and other accessories and components for use in percutaneous placement of the gastrostomy catheter.

The Cope Gastrointestinal Suture Anchor Set consists of an introducer needle with a preloaded suture anchor, a special .018/.035 inch wire guide with a spring coil tip, and a second suture anchor in a loading cannula.

### INTENDED USE

The Russell Gastrostomy Tray is intended for percutaneous placement of a balloon catheter for gastrostomy feeding.

The Cope Gastrointestinal Suture Anchor Set is intended for anchoring the anterior wall of the stomach to the abdominal wall prior to introduction of interventional catheters.

### CONTRAINDICATIONS

None known

### WARNINGS

None known

### PRECAUTIONS

- This product is intended for use by physicians trained and experienced in percutaneous placement of gastrostomy catheters. Standard techniques for placement of wire guides, introducers, and catheters should be employed.
- Do not use air to inflate the balloon; the balloon will not remain inflated when air is used.
- Manipulation of the Cope Gastrointestinal Suture Anchor Set requires fluoroscopic control.
- During placement of suture anchors, withdrawal or manipulation of the distal spring coil portion of the mandril wire guide through the needle tip may result in breakage.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### INSTRUCTIONS FOR USE

1. After administering light sedation, insert a flexible fiberoptic gastroscope into the stomach and perform a preliminary gastric examination. Direct the gastroscope tip anteriorly; dim the room lights to allow visualization of the light transmitted through the gastric and abdominal walls. (To maintain gastric distension throughout the procedure, repeatedly insufflate air.)
2. If suture anchors are to be used, place the anchors as directed in the following section.
3. After preparing the access site with antiseptic and local anesthetic, advance the introducer needle through the anterior abdominal wall and into the stomach under endoscopic visualization.
4. Introduce the wire guide through the needle and advance it into the stomach, then remove the needle, leaving the wire guide in place.
5. Adjacent to the wire guide, make a small incision through the muscle and fascia of the upper left quadrant. **NOTE:** Penetration of the peritoneum is not necessary. **NOTE:** The incision should be at least as wide as the dilator and sheath to allow their easy passage.
6. Dilate the tract through the abdominal wall and into the stomach, using the supplied dilators in sequence. Endoscopy should be used to visualize passage of the dilators into the stomach.
7. Lubricate the lumen of the supplied Peel-Away sheath, then pass the dilator/sheath assembly over the wire guide as a unit.
8. Advance the wire guide and dilator/sheath assembly into the stomach as a unit. **NOTE:** Prior incision of the skin, muscle, and fascia, as well as gastric distension via the endoscope, facilitates this procedure. In addition, counterpressure by the endoscopist or clockwise-counterclockwise rotary motion as the sheath is advanced may assist. **NOTE:** Before proceeding, both the dilator and Peel-Away sheath set must be endoscopically observed.
9. Leaving the Peel-Away sheath in the stomach, remove the dilator and wire guide.
10. Advance the lubricated gastrostomy catheter through the lubricated sheath until it is observed in the stomach.
11. Inflate the balloon with 5 ml of sterile water. Check to ensure that the

valve is functioning properly and that the balloon remains inflated.

12. **Close the red clamp on the inflation sidearm to help prevent inadvertent balloon deflation.**
13. After ensuring the system is leakproof, peel away the sheath and withdraw the balloon catheter to effect apposition of the stomach to the anterior abdominal wall.
14. Place slight tension on the catheter to create fixation of the stomach to the abdominal wall; suture at skin level.
15. Place the gastrostomy catheter on dependent drainage overnight. Begin feedings the following day.

### **Cope Gastrointestinal Suture Anchor Set**

1. Insufflate stomach with approximately 500 cc of air through a nasogastric tube.
2. Localize the puncture site and administer local anesthetic to the peritoneum. **NOTE:** The preferred puncture site into the gastric lumen is at the junction of the left anterior rib cage and the lateral border of the rectus abdominis muscle, provided there is no interposition of the liver or colon.
3. Attach a 10 ml syringe, half filled with contrast medium, to the introducer needle that has been preloaded with the suture anchor.
4. Advance the introducer needle through the anterior abdominal wall and then rapidly thrust the needle another 2-5 cm into the air-filled stomach. **Air should be freely aspirated. Injection of contrast medium under fluoroscopy will aid in confirming proper intragastric needle position.**
5. Remove the 10 ml syringe from the needle hub. While maintaining slight tension on the trailing suture, introduce the distal spring coil .035 inch portion of the wire guide into the needle and use it to push the suture anchor out of the needle into the stomach cavity. **Maintaining slight tension on the suture during anchor introduction will prevent it from bunching and lodging within the needle lumen.**
6. Remove the introducer needle over the wire guide and trailing suture needle. **The .018 inch part of the wire guide must be positioned adjacent to the suture needle to allow clearance and removal of the introducer needle.**
7. With the wire guide still in position, apply traction to the suture to pull the anterior wall of the stomach against the abdominal wall.
8. While maintaining traction on the suture anchor, secure the suture to the skin under slight tension.
9. Remove the wire guide.
10. To place the second suture anchor, reintroduce the introducer needle as described in step 4 of this section.
11. Advance the loading cannula containing the second suture anchor into the hub of the introducer needle.
12. Advance the wire guide through the cannula into the introducer needle, pushing the second suture anchor into the abdominal cavity.
13. Follow steps 6 through 8, of this section with the second suture anchor.
14. The wire guide may now be used to facilitate passage of fascial dilators and introduction of the gastrostomy catheter. Take care not to advance the entire distal spring coil segment of the wire guide into the abdominal cavity.

**NOTE:** The suture may be left in place for two weeks while tract formation occurs, or it may be cut after gastrostomy catheter placement. This releases the anchor into the stomach, allowing its passage via the gastrointestinal system.

### **HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, cool, dry place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

### **REFERENCES**

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.