- EN**  
3 **Bush SL™ and Bush DL™ Ureteral Illuminating Catheters**  
Instructions for Use
- CS**  
4 **Ureterální prosvětlovací katetry Bush SL™ a Bush DL™**  
Návod k použití
- DA**  
6 **Bush SL™ og Bush DL™ katetre til belysning af ureter**  
Brugsanvisning
- DE**  
8 **Bush SL™ und Bush DL™ Ureter-Katheter mit Beleuchtung**  
Gebrauchsanweisung
- EL**  
9 **Φωτιζόμενοι ουρητηρικοί καθετήρες Bush SL™ και Bush DL™**  
Οδηγίες χρήσης
- ES**  
11 **Catéteres de transiluminación ureteral Bush SL™ y Bush DL™**  
Instrucciones de uso
- FR**  
13 **Sondes urétérales pour transillumination Bush SL™ et Bush DL™**  
Mode d'emploi
- HU**  
15 **Bush SL™ és Bush DL™ világító ureterális katéterek**  
Használati utasítás
- IT**  
16 **Cateteri ureterali transilluminanti Bush SL™ e Bush DL™**  
Istruzioni per l'uso
- NL**  
18 **Bush SL™ en Bush DL™ ureterale illuminatiekatheters**  
Gebruiksaanwijzing
- NO**  
20 **Bush SL™ og Bush DL™ gjennomlysningskatetre for ureter**  
Bruksanvisning
- PL**  
22 **Cewniki moczowodowe Bush SL™ i Bush DL™ do transiluminacji**  
Instrukcja użycia
- PT**  
23 **Cateteres para transiluminação ureteral Bush SL™ e Bush DL™**  
Instruções de utilização
- SV**  
25 **Bush SL™ och Bush DL™ uretärkatetrar med belysning**  
Bruksanvisning

## BUSH SL™ and BUSH DL™ URETERAL ILLUMINATING CATHETERS

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

### DEVICE DESCRIPTION

#### **Bush SL sets contain:**

- Two 5.0 French single-lumen transparent vinyl catheters, 80 cm long, with radiopaque and non-radiopaque markings
- Two fiberoptic light guides with anodized aluminum plugs

#### **Bush DL sets contain:**

- Two 7.0 French dual-lumen transparent vinyl catheters, 80 cm long, with radiopaque and non-radiopaque markings
- Two fiberoptic light guides with anodized aluminum plugs
- Two 7.0 French polyvinylchloride drainage tubes, 29 cm long
- Two 10.0 French polyvinylchloride connecting tubes, 30 cm long

**NOTE:** Set components may vary.

### PRODUCT RECOMMENDATIONS

The anodized aluminum plugs on the light guides fit a standard CIRCON-A.C.M.I. light source. If using a different light source, an adapter may be required to successfully use the Bush DL or Bush SL Ureteral Illuminating Catheter. Adapters are available from various light source manufacturers.

### INTENDED USE

Used for transillumination of the ureters during laparoscopic or open surgical procedures. Transillumination helps to identify, and minimize potential for trauma to, the ureters.

### CONTRAINDICATIONS

None known

### WARNINGS

None known

### PRECAUTIONS

- The anodized aluminum plug conducts heat. Allow the plug to cool down prior to attempting to unplug.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### INSTRUCTIONS FOR USE:

1. Before beginning procedure, ensure that the light source you are using accepts the anodized aluminum plugs included with this product, or that you have an appropriate adapter available.
2. Remove the Fitting from the proximal end of the catheter.
3. Place the catheter cystoscopically. On standard ureteral lengths of 24-26 cm, position the 25 cm ink mark (second from distal tip) at the ureteral vesicle junction (UVJ). The catheter position may need to be adjusted on longer or shorter length ureters. **NOTE:** Radiopaque bands permit cystoscopic and fluoroscopic confirmation of catheter placement, as well as the location of the illuminated area. Improper positioning of the catheter may impair illumination of the ureter.
4. Disconnect the light fiber extension and remove the scope from the catheter.

5. If illuminating both ureters, repeat steps 1-3 on the opposite ureter. **NOTE:** The catheters provided in this set are marked with differently colored bands (red and black) to differentiate between the right and left ureters.
6. Reinsert the Fitting into the proximal end of the catheter.
7. Insert the light plug into the fiberoptic light source.  
**NOTE:** A poor or improper light source connection may result in a decreased amount of light emanating from the catheter.  
**NOTE:** This catheter will not transmit thermal energy along its light fibers to patient tissue. Any excessive thermal damage from the light source will manifest at the catheter/light plug junction, and will not be transmitted along the catheter length.  
**NOTE:** If using Bush DL catheters, the drainage tubes may be connected to a collection bag at this point if necessary.
8. Activate the light source to transilluminate the ureter.  
**NOTE:** Start illumination with the light source at the lowest setting, as many light sources produce thermal energy at varying temperatures. This will limit the possibility of thermal damage at the catheter/light plug junction. It is inadvisable to use any light source at its highest setting, unless the light source's actual thermal energy output is known.  
**NOTE:** Due to anatomical differences, visualization of transilluminated ureters may vary between patients. Placing the laparoscope in a lateral port site may enhance visualization in many cases. For maximum transillumination of the ureter(s), reduce the intensity of the laparoscope's light or, in an open case, reduce the operative field lights.  
**NOTE:** High-energy light sources such as Xenon may cause overheating of the anodized aluminum plug. An appropriate adapter (available from most light source manufacturers) will ensure product safety and functionality. For Storz Xenon light source, use "Turret Adapter 487UO."

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.