EN 3	Foley Catheter Tip Punch Instructions for Use
CS 5	Děrovačka hrotu pro katetry Foley Návod k použití
DA 7	Stanse til Foley kateterspidser Brugsanvisning
DE 9	Stanze für Foley-Katheter-Spitzen Gebrauchsanweisung
EL 12	Διατρητής άκρου καθετήρα Foley Οδηγίες χρήσης
ES 14	Sacabocados para puntas de catéteres Foley Instrucciones de uso
FR 17	Perforateur d'extrémité pour sonde de Foley Mode d'emploi
HU 19	Foley katétervég-lyukasztó Használati utasítás
IT 21	Punzone per punta di catetere Foley Istruzioni per l'uso
NL 24	Pons voor Foley-kathetertip Gebruiksaanwijzing
NO 26	Hullemaskin for Foley kateterspisser Bruksanvisning
PL 29	Przebijak końcówki cewnika Foley'a Instrukcja użycia
PT 31	Punção para ponta de cateter Foley Instruções de utilização
SV 33	Stans för Foley-kateterspets Bruksanvisning

FOLEY CATHETER TIP PUNCH

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

INTENDED USE

Used for placing an end hole in Foley catheters sized 14 French or larger.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

None known

INSTRUCTIONS FOR USE

1. Place Foley catheter tip on stem of punch.
2. Screw plunger completely through tip.
3. Remove Foley catheter and discard tip material.

Disassembly and Reassembly

To disassemble the Foley Catheter Tip Punch, unscrew the thumb screw by twisting it counterclockwise until the thumb screw comes completely out of the body of the tip punch fixture. Do not attempt to remove the punch from the fixture.

To reassemble the Foley Catheter Tip Punch, carefully reinsert the thumb screw by twisting it clockwise until the thumb screw sits securely within the fixture. Do not over tighten the thumb screw, as this can damage the device.

Cleaning, Disinfection and Sterilization

NOTE: All cleaning, disinfection and sterilization methods listed herein have been validated.

To Clean

1. Disassemble the Foley Catheter Tip Punch per the instructions listed within the IFU.
2. Remove all residual organic matter, blood and irrigation solution to allow the sterilization or disinfection medium to contact the surface of the device. **NOTE:** Ethylene oxide gas and disinfection solutions do not effectively penetrate residual organic matter.
3. Prepare a mild, neutral-pH enzymatic detergent (i.e., Enzol using 1 ounce per gallon) and warm (30-43°C) tap or Deionized (DI) water. Immerse disassembled device and soak for 20 minutes. Use a soft brush to remove all residual organic matter, blood and irrigation solution from the device.
4. Prepare a fresh, neutral-pH enzymatic detergent (i.e., Enzol using 1 ounce per gallon) and warm (30-43°C) tap or DI water. Immerse disassembled device and perform ultrasonic cleaning for 20 minutes.
5. Rinse thoroughly with clean, running tap or DI water.
6. Dry all exterior surfaces with a soft cloth.
7. Dry all openings with an air pistol (compressed air).
8. Prior to disinfection/sterilization and reuse, this product should be inspected and reassembled to verify that it functions properly. Grasping edges and all surfaces should be smooth and free of abrasions. The device must be disassembled again prior to disinfection and/or sterilization.

NOTE: Automated cleaning methods are not compatible with the Foley Catheter Tip Punch.

To Disinfect

NOTE: CIDEX® Solution is a high-level disinfection solution that has been found to be compatible with the Foley Catheter Tip Punch.

1. Disassemble the Foley Catheter Tip Punch per the instructions listed within the IFU, if the device has been reassembled between cleaning and disinfection.
2. Activate CIDEX® Solution per the manufacturer's instructions. Test the activated solution prior to use with CIDEX® Solution Test Strips.
3. Immerse disassembled device completely, filling all hard-to-reach areas and eliminating air pockets, in CIDEX® Activated Dialdehyde Solution for a minimum of 45 minutes at 25°C.
4. Rinse thoroughly in a room temperature Sterile Water for Injection (SWFI) water bath. Rinse by immersing device completely for a minimum of 1 minute. Manually flush all hard-to-reach areas with SWFI water.
5. Agitate device under SWFI water, bring above water level, then re-immerses.
6. Repeat steps 4 and 5, 2 additional times using fresh SWFI for each of the 3 rinses.

NOTE: Automated disinfection methods are not compatible with the Foley Catheter Tip Punch.

To Sterilize

NOTE: Autoclave sterilization is not recommended for the Foley Catheter Tip Punch.

This device can be sterilized by ethylene oxide (EO) gas. EO sterilization should be performed in accordance with industry standards with a validated sterilization cycle. The device must be disassembled prior to sterilization. The recommended cycle for sterilizing this device using 100% EO in a double-pouched configuration is:

Temperature = 54.4°C

Humidity = 50-80%RH

EO concentration = 735 mg/L EO

Gas dwell time = minimum 60 minutes

LIMITED WARRANTY STATEMENT

Cook Incorporated ("Company") warrants to the purchaser of this product (the "Product") that at the time of manufacture, the Product was manufactured and tested in accordance with Good Manufacturing Practices specified by the United States Food and Drug Administration and was true to label. In the unlikely event of a Product failure due to defects in material or workmanship within the warranty period, the Product will be repaired or replaced, at Cook's option, free of charge.

Because the Company has no control over the conditions under which this Product is used, diagnosis of the patient, anatomical or biological differences among patients, method of use or administration, or handling of the Product after it leaves the Company's possession, the Company does not warrant a good effect or against an ill effect following the Product's use. The foregoing warranty is exclusive and in lieu of all other warranties whether written, oral, or implied. There are no warranties of merchantability or fitness for a particular purpose. No representative of the Company may change any of the foregoing, and the purchaser accepts the Product subject to all terms hereof. The Company reserves the right to change or discontinue this Product without notice.

Limits and Exclusions

The warranty does not apply in the event of misuse, negligence, or accidental damage of Product. Warranty services will be provided only to the original purchaser at the location where the Product was originally installed. The warranty may not be transferred without the prior written consent of the Company.

Warranty Period

This limited warranty begins on the date the Product ships from the Company's facilities and ends 6 months later.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Reuseable. Sterile if the package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.