EN
3

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T _ J - B G C _ R E V O

WORD CATHETER – SILICONE BARTHOLIN GLAND BALLOON SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- *Set components may vary*
- Silicone Balloon Catheter
- Syringe
- Scalpel

INTENDED USE

Used for treatment of abscesses and cysts of the Bartholin gland.

CONTRAINDICATIONS

None Known

WARNINGS

- Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas.
- Do not overinflate. Using excessive pressure to inflate the balloon on this device can cause the balloon to rupture.
- Consider underlying malignancy in women over age 40 who present with Bartholin Gland cysts or abscesses.

PRECAUTIONS

- Do not exceed the maximum 3ml inflation volume for this balloon device.
- Do not use petroleum based ointments or lubricants since they may damage silicone and burst balloon.
- Epithelization time may vary depending on each patient. The Word Catheter-Silicone Bartholin Gland Balloon is not intended to be left indwelling for periods of time greater than twenty-eight (28) days.

INSTRUCTIONS FOR USE

Word Catheter – Silicone Bartholin Gland Balloon Set:

1. Clean area surrounding Bartholin gland with appropriate antiseptic solution.
2. Using enclosed scalpel, make an incision in the outer wall of the cyst, preferably inside the hymenal ring, to break up any loculations.
3. Once the abscess has drained, insert the deflated silicone catheter through the incision and into the cystic body.
4. Using the enclosed syringe, carefully inflate the silicone catheter with sterile saline until the balloon is sufficiently retained in the cyst, not exceeding maximum inflation volume of the device (3ml).
5. Once inflated, remove the syringe and needle, leaving the saline-inflated Word catheter balloon within the cystic mass.
6. Tuck free end of the catheter up into the vagina.
7. After epithelization of a new orifice is accomplished, carefully deflate the balloon using a sterile syringe with needle tip, and remove the catheter.

NOTE: Epithelization time may vary depending on each patient. The Word Bartholin Gland Catheter is not intended to be left indwelling for periods of time greater than twenty-eight (28) days.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

Placement of a Word Catheter: A resident training model. Owen J, Koza J, Shiblee T, Robertazzi R, Cabbad M, Hsu C. American Journal of Obstetrics and Gynecology. Volume 192, Issue 5, May 2005.

Office Management of Bartholin Gland Cysts and Abscesses. Hill DA, Lense JJ. American Family Physician. Volume 57, Issue 7, April 1, 1998.

Management of Bartholin Duct Cysts and Abscesses, A systematic Review. Wechter ME, Wu JM, Marzano D, Haefner H. Obstetrical and Gynecological Survey, Volume 64, Number 6, June 2009.

Draining a cyst or abscess in a Bartholin gland with a Word Catheter. Scott P. JAAPA; Official Journal of the American Academy of Physician Assistants. Volume 16, Issue 12, 2003.

The simple outpatient management of Bartholin's Abscess using the Word catheter: A preliminary study. Haider Z, Condous G, Kirk E, Mukri F, Bourne T. The Australian and New Zealand Journal of Obstetrics and Gynecology. Volume 47, Issue 2, 2007.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.