EN  
3**Laparoscopic and Hysteroscopic Myoma Screws**

Instructions for Use

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Bruksanvisning

T \_ J - L H M S \_ R E V 1

## LAPAROSCOPIC AND HYSTEROSCOPIC MYOMA SCREWS

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### INTENDED USE

Used to stabilize and manipulate non-degenerating myomas or fibroids, or other tissues during laparoscopic and hysteroscopic procedures.

### CONTRAINDICATIONS

Use of these devices is contraindicated in the presence of infection, pregnancy or a degenerating myoma or fibroid.

### WARNINGS

None known

### PRECAUTIONS

- The Laparoscopic Myoma Screw and the Hysteroscopic Myoma Screw are intended for use by physicians trained and experienced in laparoscopic and hysteroscopic procedures.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### INSTRUCTIONS FOR USE

#### Laparoscopic Myoma Screw

1. Insert the laparoscopic myoma screw into the abdomen through a 5 mm or larger access port.
2. Place distal tip in contact with tissue to be anchored.
3. Anchor device by turning the handle clockwise and applying slight pressure.
4. Manipulate and/or dissect attached tissue.
5. Remove device from patient and dispose of properly.

#### Hysteroscopic Myoma Screw

**CAUTION: Deflection of hysteroscope can damage this device.**

1. Insert device into a 2 mm or larger scope channel.
2. Place distal tip in contact with tissue to be anchored.
3. Anchor device by turning the handle clockwise and applying slight pressure.
4. Manipulate and/or dissect attached tissue.
5. Remove device from scope and dispose of properly.

### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.