



Margolin

HSG CANNULA

Achieve enhanced control during placement with removable stylet.

Easily placed and especially effective for a stenotic cervix.

Balloonless design enables visualization of the uterus with no artifacts.

Soft silicone wedge is atraumatic.





Margolin

HSG CANNULA



Equipped with a 9.0 French tip on a soft silicone wedge, the Margolin is helpful in accessing the uterine cavity for patients with a stenotic cervical canal. It is easily placed and features a balloonless design that enables visualization of the uterus without potential artifacts.

Order Number	Reference Part Number	Fr	Length cm
G16504	J-MHSG-442500	16	25



Be the first to learn about Cook Medical's new products, upcoming events and live webcasts. Sign up for e-newsletters at <http://www.cookmedical.com/wh/products.do>.

Other Imaging Catheters



Silicone Balloon HSG Catheter

Used to perform both SHG and HSG, this catheter has a small diameter that allows easier passage through the cervix and a small balloon that enables visualization of the uterus and fallopian tubes. Sold in boxes of 10.

Order Number	Reference Part Number	Fr	Length cm	Balloon Volume mL	Syringe Volume mL
Closed Tip with Sideport					
G17472	J-CHSG-503001	5	30	1	1
G17471	J-CHSG-703001	7	30	1.5	3
Open Tip					
G17147	J-CHSG-503000	5	30	1	1
G17870	J-CHSG-703000	7	30	1.5	3



Sholkoff Balloon Hysterosalpingography Catheter

Used to inject contrast medium into the uterine cavity and fallopian tube(s) for HSG or evaluation of tubal patency, the Sholkoff features a radiopaque introducer sheath that facilitates placement of the catheter into the uterine cavity.

Order Number	Reference Part Number	Fr	Length cm	Balloon Volume mL	Syringe Volume mL
G16575	J-SBH-683000	6.8	30	1.4	3

IN THIS BROCHURE

Products shown intended for one-time use.

Not all part numbers shown may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

- AORTIC INTERVENTION
- CRITICAL CARE
- ENDOSCOPY
- INTERVENTIONAL RADIOLOGY
- LEAD MANAGEMENT
- PERIPHERAL INTERVENTION
- SURGERY
- UROLOGY
- WOMEN'S HEALTH

EMEA: EDI – www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland – French: +41 448009609, fr.orders@cookmedical.com
 Switzerland – Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland – German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI – www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
 E-mail: cau.custserv@cookmedical.com

